

ATAD1 antibody

Catalog Number: 108276

Product name

ATAD1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ATAD1 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in rat brain tissue. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of ATAD1 recombinant protein (Accession Number: BC010868). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108276(ATAD1 Antibody) at dilution of 1:300



IP Result of anti-ATAD1 (IP:Catalog No:108276, 3ug; Detection:Catalog No:108276 1:300) with mouse brain tissue lysate 4000ug.