

# STAR antibody

Catalog Number: 115676

## Product name

STAR antibody

## Specificity

Human, Mouse, Rat, Pig; other species not tested.

## Antibody description

STAR Rabbit Polyclonal antibody. Positive WB detected in human adrenal gland tissue, pig adrenal gland tissue. Positive IHC detected in human ovary tissue, mouse testis tissue. Observed molecular weight by Western-blot: 32 kDa

## Preparation

This antibody was obtained by immunization of STAR recombinant protein (Accession Number: NM\_000349). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



human adrenal gland tissue were subjected to SDS PAGE followed by western blot with Catalog No:115676(STAR antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded human ovary tissue slide using Catalog No:115676(STAR Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human ovary tissue slide using Catalog No:115676(STAR Antibody) at dilution of 1:50

# STAR antibody



Catalog Number: 115676

---

(under 40x lens)