

UCK2 antibody

Catalog Number: 116676

Product name

UCK2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

UCK2 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 33-35 kDa

Preparation

This antibody was obtained by immunization of UCK2 recombinant protein (Accession Number: NM_012474). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116676(UCK2 antibody) at dilution of 1:500



IP Result of anti-UCK2 (IP:Catalog No:116676, 4ug; Detection:Catalog No:116676 1:500) with Jurkat cells lysate 4000ug.