

ESAM antibody



Catalog Number: 110376

Product name

ESAM antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

ESAM Rabbit Polyclonal antibody. Positive WB detected in human lung tissue, human heart tissue. Observed molecular weight by Western-blot: 56-58kd

Preparation

This antibody was obtained by immunization of ESAM recombinant protein (Accession Number: NM_138961). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:110376(ESAM antibody) at dilution of 1:500