

Collagen Type II antibody

Catalog Number: 109376

Product name

Collagen Type II antibody

Specificity

Human; other species not tested.

Antibody description

Collagen Type II Rabbit Polyclonal antibody.
Positive WB detected in transfected cells. Positive IHC detected in human lung tissue.

Preparation

This antibody was obtained by immunization of Collagen Type II recombinant protein (Accession Number: BC007252). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:250-1:1000

IHC: 1:20-1:200

Validations




WB result of anti-COL2A1 (Catalog No:109376) with overexpressed protein COL2A1(1-268aa)-



GFP.

Immunohistochemistry of paraffin-embedded human lung tissue slide using Catalog No:109376(COL2A1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human lung tissue slide using Catalog No:109376(COL2A1 Antibody) at dilution of 1:50 (under 40x lens)