

Human HDAC4 (aa 612-1084) recombinant protein



Catalog Number: 504776

General Information

Protein Construction

A DNA sequence encoding the human HDAC4 (Met612-Leu1084) was expressed and purified with two additional amino acids (Gly & Pro) at the N-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The secreted recombinant human HDAC4 consists of 475 amino acids and predicts a molecular mass of 50.9 KDa. The apparent molecular mass of the

protein is approximately 51 KDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile 20mM Tris,500mM NaCl, pH 7.4, 10%gly.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human HDAC4 Protein (aa 612-1084) SDS-PAGE