

DDR1 antibody

Catalog Number: 109776

Product name

DDR1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DDR1 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue, human kidney tissue.

Preparation

This antibody was obtained by immunization of DDR1 recombinant protein (Accession Number: NM_001387917). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC


Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations

Immunohistochemical of paraffin-embedded human testis using Catalog No:109776(DDR1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human testis using Catalog No:109776(DDR1 antibody) at dilution of 1:50 (under 40x lens)