

# Anti-LAIR1 antibody

Catalog Number: 141086

## Product name

Anti-LAIR1 antibody

## Specificity

Human LAIR1/LAIR-1/CD305

## Antibody description

Rabbit Polyclonal to Human LAIR1

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Human LAIR1/LAIR-1/CD305, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 5-20 µg/ml

## Validations



Items	Lanes	A
Sample		mouse liver tissue lysate
Sample Volume (µg/lane)		30
Gel		13% SDS-PAGE reducing gel
Recommended Concentration		5.20 µg/ml
Secondary Antibody		Dylight 800-labeled Antibody to Rabbit IgG (H+L) at 1:5000 dilution.
Developed using Odyssey imaging system.		
Explanation		Predicted band size : 31 kDa Observed band size : 30 kDa

Human LAIR1/LAIR-1/CD305 Western blot (WB)  
15857