

GLYCTK antibody



Catalog Number: 111086

Product name

GLYCTK antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

GLYCTK Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer. Positive WB detected in human lung tissue, human liver tissue. Observed molecular weight by Western-blot: 55kd

Preparation

This antibody was obtained by immunization of GLYCTK recombinant protein (Accession Number: BC021896). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

Validations



human lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:111086(GLYCTK antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:111086(GLYCTK antibody) at dilution of 1:50 (under 10x lens)