

NDUFV1 antibody



Catalog Number: 113086

Product name

NDUFV1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NDUFV1 Rabbit Polyclonal antibody. Positive IP detected in A431 cells. Positive WB detected in A431 cells, human brain tissue. Positive IHC detected in human lung cancer tissue. Observed molecular weight by Western-blot: 51 kDa

Preparation

This antibody was obtained by immunization of NDUFV1 recombinant protein (Accession Number: NM_007103). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

IP: 1:200-1:2000

IHC: 1:20-1:200

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:113086(NDUFV1 antibody) at dilution of 1:500

Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:113086(NDUFV1 antibody) at dilution of 1:50 (under 10x lens)

IP Result of anti-NDUFV1 (IP:Catalog No:113086, 4ug; Detection:Catalog No:113086 1:500) with A431 cells lysate 2000ug.