

PPAT antibody

Catalog Number: 114086

Product name

PPAT antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PPAT Rabbit Polyclonal antibody. Positive IP detected in HepG2 cells. Positive WB detected in human stomach tissue, HEK-293 cells, HepG2 cells. Observed molecular weight by Western-blot: 50-56kd

Preparation

This antibody was obtained by immunization of PPAT recombinant protein (Accession Number: NM_002703). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



human stomach tissue were subjected to SDS PAGE followed by western blot with Catalog No:114086(PPAT antibody) at dilution of 1:500



IP Result of anti-PPAT (IP:Catalog No:114086, 4ug; Detection:Catalog No:114086 1:500) with HepG2 cells lysate 2400ug.