

MOGS Polyclonal Antibody



Catalog Number: 166086

Product name

MOGS Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to MOGS

Preparation

Antigen: A synthetic Peptide of human MOGS

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB

Dilutions

WB 1:500 - 1:1000

Validations

Mouse liver



Western blot - MOGS Polyclonal Antibody

Western blot analysis of extracts of mouse liver, using MOGS antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.