

TMED3 antibody



Catalog Number: 116086

Product name

TMED3 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

TMED3 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, HeLa cells, MCF7 cells, mouse skin tissue, mouse small intestine tissue, HT-1080 cells. Observed molecular weight by Western-blot: 21-25kd

Preparation

This antibody was obtained by immunization of TMED3 recombinant protein (Accession Number: NM_007364). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116086(TMED3 antibody) at dilution of 1:1000