

# Anti-eIF4E (8C8) Mouse antibody



Catalog Number: 168086

## Product name

Anti-eIF4E (8C8) Mouse antibody

## Specificity

Human, Mouse, Rat

## Antibody description

eIF4E (8C8) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human eif4e protein expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

WB


## Dilutions

WB: 1/1000

## Validations



Western blot detection of eIF4E in Jurkat, HeLa, MCF7, K562, C6, 3T3, PC12, HepG2, CHO-K1, Raji, Raw264.7 and Ramos cell lysates using eIF4E mouse mAb (1:2000 diluted). Predicted band size: 25KDa. Observed band size: 25KDa.



Western blot detection of eIF4E in MDA-MB-468, ZR75-1, MCF7, T47D and CHO-K1 cell lysates using eIF4E mouse mAb (1:1000 diluted). Predicted band size: 25KDa. Observed band size: 25KDa.