

Anti-eIF4E (8C8) Mouse antibody



Catalog Number: 168086

Product name

Anti-eIF4E (8C8) Mouse antibody

Specificity

Human, Mouse, Rat

Antibody description

eIF4E (8C8) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human eIF4E protein expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB


Dilutions

WB: 1/1000

Validations



Western blot detection of eIF4E in Jurkat, HeLa, MCF7, K562, C6, 3T3, PC12, HepG2, CHO-K1, Raji, Raw264.7 and Ramos cell lysates using eIF4E mouse mAb (1:2000 diluted). Predicted band size: 25kDa. Observed band size: 25kDa.



Western blot detection of eIF4E in MDA-MB-468, ZR75-1, MCF7, T47D and CHO-K1 cell lysates using eIF4E mouse mAb (1:1000 diluted). Predicted band size: 25kDa. Observed band size: 25kDa.