

# Human FGF2/FGFB recombinant protein



Catalog Number: 518086

## General Information

### Gene Name Synonym

Fibroblast growth factor 2, FGF-2, Basic fibroblast growth factor, Bfgf, Heparin-binding growth factor 2, HBGF-2, FGF2, FGFB

### Protein Construction

A DNA sequence encoding the Human FGF2/FGFB Met134-Ser288 is expressed.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL method.

### Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

## Formulation

Lyophilized from a 0.2  $\mu$ m filtered solution of 20mM TrisHCl, 150mM NaCl, pH 7.5.

## Usage Guide


### Storage

Lyophilized protein should be stored at <-20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days.

### Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 $\mu$ g/ml. Dissolve the lyophilized protein in ddH<sub>2</sub>O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## SDS-PAGE



Human FGF2/FGFB recombinant protein