

Anti-JAB1 Rabbit antibody



Catalog Number: 169086

Product name

Anti-JAB1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

JAB1 Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human JAB1

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP

Dilutions


WB: 1/500-1/2000

IHC: 1/50-1/200

ICC/IF: 1/50-1/200

IP: 1/50

Validations



Western blot analysis of JAB1 expression in MCF7 cell lysate.