

# HMGN1 antibody



Catalog Number: 111486

## Product name

HMGN1 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

HMGN1 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, HeLa cells. Positive IF detected in HepG2. Observed molecular weight by Western-blot: 17kd

## Preparation

This antibody was obtained by immunization of HMGN1 recombinant protein (Accession Number: NM\_004965). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:111486(HMGN1 antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using HMGN1 antibody Catalog No:111486 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).