

Anti-E Cadherin (7H12) Mouse antibody



Catalog Number: 168486

Product name

Anti-E Cadherin (7H12) Mouse antibody

WB: 1/500 - 1/2000

IHC: 1/200 - 1/1000

Specificity

Human, Mouse, Monkey

FC: 1/200 - 1/400

ELISA: 1/10000


Antibody description

E Cadherin (7H12) Mouse monoclonal antibody

Validations

Preparation

Antigen: Purified recombinant fragment of human CDH1 expressed in E. Coli.



Formulation

Ascitic fluid containing 0.03% sodium azide.

Western blot analysis of extract from T47D cells using E-Cadherin Mouse mAb at 1:1000 dilution. Predicted band size: 135KDa. Observed band size:

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG1

Applications

WB, IHC-P, FC

Dilutions



135KDa.

Western blot detection of E-Cadherin in C6 and HaCat cell lysates using E Cadherin mouse mAb (dilution 1:500). Predicted band size: 135KDa. Observed band size: 135KDa.