

JMJD3 antibody



Catalog Number: 111886

Product name

Polyclonal

JMJD3 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse; other species not tested.

Applications

Antibody description

ELISA, WB

JMJD3 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 180 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001080424). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111886(JMJD3 Antibody) at dilution of 1:600

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality