

PRAME Polyclonal Antibody



Catalog Number: 164886

Product name

PRAME Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to PRAME

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human PRAME (NP_996839.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:20 - 1:50

Validations



Immunofluorescence - PRAME Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using PRAME antibody . Blue: DAPI for nuclear staining.