

# Anti-ARFRP1/ARP Rabbit antibody



Catalog Number: 168886

## Product name

Rabbit IgG

Anti-ARFRP1/ARP Rabbit antibody

## Applications

## Specificity

WB, ICC/IF, IP

Human, Mouse, Rat

## Dilutions

## Antibody description

1:1000-1:5000

ARFRP1/ARP Rabbit monoclonal antibody

1:20-1:50

## Preparation

1:20-1:50

Antigen: A synthetic peptide of human ARFRP1

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot detection of ARFRP1/ARP in C6 cell lysates using ARFRP1/ARP Rabbit mAb(1:1000 diluted).Predicted band size:23kDa.Observed band size:23kDa.