

# Cathepsin A antibody



Catalog Number: 108886

---

## Product name

Cathepsin A antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Cathepsin A Rabbit Polyclonal antibody. Positive WB detected in Sp2/0 cells. Observed molecular weight by Western-blot: 54-60 kDa

## Preparation

This antibody was obtained by immunization of Cathepsin A recombinant protein (Accession Number: NM\_001127695). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



Sp2/0 cells were subjected to SDS PAGE followed by western blot with Catalog No:108886(CTSA antibody) at dilution of 1:300