

Human UBE2S/E2EPF/OK/SW-cl.73 (GST tag) recombinant protein



Catalog Number: 519886

General Information

Gene Name Synonym

Ubiquitin-Conjugating Enzyme E2 S, E2-EPF, Ubiquitin Carrier Protein S, Ubiquitin-Conjugating Enzyme E2-24 kDa, Ubiquitin-Conjugating Enzyme E2-EPF5, Ubiquitin-Protein Ligase S, UBE2S, E2EPF

Protein Construction

A DNA sequence encoding the Human UBE2S/E2EPF/OK/SW-cl.73 Met1-Leu222 is expressed with a GST tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human UBE2S/E2EPF/OK/SW-cl.73 (GST tag)
recombinant protein