

EIF2C3-Specific antibody



Catalog Number: 110186

Product name

EIF2C3-Specific antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

EIF2C3-Specific Rabbit Polyclonal antibody. Positive WB detected in SKOV-3 cells, HepG2 cells. Observed molecular weight by Western-blot: 90-100 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_024852). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



SKOV-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:110186(EIF2C3-Specific Antibody) at dilution of 1:300