

TES Polyclonal Antibody



Catalog Number: 165186

Product name

TES Polyclonal Antibody

Specificity

Human, Rat

Antibody description

Polyclonal antibody to TES

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human TES (NP_056456.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB, IHC

Dilutions


WB 1:500 - 1:2000
IHC 1:50 - 1:200

Validations



Immunohistochemistry - TES Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat ovary using TES antibody at dilution of 1:100 (40x



lens).

Immunohistochemistry - TES Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung cancer using TES antibody at dilution



of 1:100 (40x lens).


Immunohistochemistry - TES Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver cancer using TES antibody at dilution

TES Polyclonal Antibody



Catalog Number: 165186



Immunohistochemistry - TES Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human placenta using TES antibody at dilution of 1:100 (40x lens).

of 1:100 (40x lens).