

Anti-OR4C16 antibody



Catalog Number: 176186

Product name

Anti-OR4C16 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to OR4C16

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or4c16. aa range:261-310

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.; Western Blot analysis of SH-SY5Y cells using Olfactory receptor 4C16 Polyclonal Antibody diluted at 1: 500