

Icam1 antibody

Catalog Number: 111586

Product name

Icam1 antibody

Specificity

Mouse,Rat; other species not tested.

Antibody description

Icam1 Rabbit Polyclonal antibody. Positive IHC detected in mouse liver tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 85-95 kDa

Preparation

This antibody was obtained by immunization of Recombinant Protein (Accession Number: NM_010493). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

IHC: 1:20-1:200


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111586(Icam1 antibody) at dilution of 1:4000



Immunohistochemistry of paraffin-embedded mouse liver tissue slide using Catalog No:111586(Icam1 Antibody) at dilution of 1:200 (under 10x lens)



Immunohistochemistry of paraffin-embedded mouse liver tissue slide using Catalog No:111586(Icam1 Antibody) at dilution of 1:200 (under 40x lens)