

# IFT57 Polyclonal Antibody



Catalog Number: 166586

## Product name

IFT57 Polyclonal Antibody

WB 1:500 - 1:2000

IHC 1:50 - 1:200

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to IFT57

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 170-429 of human IFT57 (NP\_060480.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications

WB, IHC


## Dilutions

## Validations



Immunohistochemistry - IFT57 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human uterine cancer using IFT57 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - IFT57 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat ovary using IFT57 antibody at dilution of 1:100 (40x lens).