

CRYGC antibody



Catalog Number: 109586

Product name

CRYGC antibody

Clonality

Polyclonal

Immunogen

[Human CRYGC Recombinant protein \(GST tag\)](#)

Ig Type

Rabbit IgG

Specificity

human, mouse, rat; other species not tested.

Applications

ELISA, WB

Antibody description

CRYGC Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 21kd

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of CRYGC recombinant protein (Accession Number: NM_020989). Purification method: Antigen affinity purified.


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109586(CRYGC antibody) at dilution of 1:500