

# KIAA0895 antibody



Catalog Number: 111986

## Product name

KIAA0895 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

KIAA0895 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, HEK-293 cells. Observed molecular weight by Western-blot: 45, 60 kDa

## Preparation

This antibody was obtained by immunization of KIAA0895 recombinant protein (Accession Number: BC110422). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:111986(KIAA0895 Antibody) at dilution of 1:300