

TAF3 antibody

Catalog Number: 115986

Product name

TAF3 antibody

Specificity

Human; other species not tested.

Antibody description

TAF3 Rabbit Polyclonal antibody. Positive WB detected in cells. Observed molecular weight by Western-blot: 140 kDa

Preparation

This antibody was obtained by immunization of TAF3 recombinant protein (Accession Number: BC024661). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: N/A

Validations



WB results of TAF3 (Catalog No:115986) with siRNA sample. (140 kDa band absent in siRNA sample)