

XIAP antibody

Catalog Number: 116986

Product name

XIAP antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

Human, Rat; other species not tested.

Antibody description

XIAP Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, C6 cells, COLO 320 cells, Jurkat cells. Positive IF detected in MCF-7 cells. Positive IHC detected in human pancreas cancer tissue. Observed molecular weight by Western-blot: 45-57 kDa

Preparation

This antibody was obtained by immunization of XIAP recombinant protein (Accession Number: NM_001204401). Purification method: Protein A purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116986(XIAP antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:116986(XIAP antibody) at dilution of 1:200 (under 40x lens)



Immunofluorescent analysis of MCF-7 cells, using

XIAP antibody



Catalog Number: 116986

XIAP antibody Catalog No: at 1:25 dilution and

Rhodamine-labeled goat anti-rabbit IgG (red).