

# Anti-NOP10 Rabbit antibody



Catalog Number: 167986

## Product name

Anti-NOP10 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

NOP10 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human NOP10

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF


## Dilutions

WB: 1/1000

IHC: 1/50

ICC/IF: 1/50

## Validations



Western blot detection of NOP10 in HeLa,HL-60  
using NOP10 Rabbit mAb(1:1000 diluted)