

Anti-ILF3/NF90 Rabbit antibody



Catalog Number: 168986

Product name

Anti-ILF3/NF90 Rabbit antibody

Specificity

Human, Mouse

Antibody description

ILF3/NF90 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human ILF3

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/2000-1/20000


IHC: 1/50

ICC/IF: 1/100

FC: 1/20

IP: 1/20

Validations



Western blot detection of ILF3/NF90 in Jurkat, C6, HeLa cell lysates using ILF3/NF90 Rabbit mAb(1:1000 diluted). Predicted band size:95kDa. Observed band size:95kDa.