

# TNFAIP3 antibody



Catalog Number: 116286

## Product name

TNFAIP3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TNFAIP3 Rabbit Polyclonal antibody. Positive WB detected in RAW 264.7 cells, HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 90 kDa, 37-40 kDa

## Preparation

This antibody was obtained by immunization of TNFAIP3 recombinant protein (Accession Number: NM\_001270507). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



RAW264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:116286(TNFAIP3 antibody) at dilution of 1:1000