

Human FLT3L / Flt3 ligand (His Tag) recombinant protein



Catalog Number: 500686

General Information

Gene Name Synonym

SL cytokine

Protein Construction

A DNA sequence encoding the human FLT3LG (P49771-1) (Thr27-Pro185) was expressed, with a polyhistidine tag at the N-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant human human FLT3LG consists

of 175 amino acids and predicts a molecular mass of 20.2 KDa. It migrates as an approximately 27 KDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 20mM Tris, 500mM NaCl, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human FLT3L / Flt3 ligand / FLT3LG Protein (His Tag) SDS-PAGE