

Mouse TIGIT/VSTM3 (His Tag) recombinant protein



Catalog Number: 501686

General Information

Protein Construction

A DNA sequence encoding the extracellular domain of mouse TIGIT (NP_001139797.1) (Met 1-Gly 141) was expressed, with a C-terminal polyhistidine tag.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Activity

Immobilized mouse TIGIT-His at 10 µg/ml (100 µl/well) can bind mouse PVR-Fch (Cat:504218), The EC₅₀ of mouse PVR-Fch (Cat:504218) is 0.31-0.73 µg/ml.

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gly 26

Molecular Mass

The secreted recombinant mouse TIGIT comprises 127 amino acids and has a calculated molecular mass of 14.2 kDa. As a result of glycosylation, the apparent molecular mass of rmSIRPA is approximately 20-25 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse TIGIT / VSTM3 Protein (His Tag) SDS-PAGE