

# Anti-SIGLEC10 antibody



Catalog Number: 175686

## Product name

Anti-SIGLEC10 antibody

## Specificity

Human, Mouse, Rat, Dog, Pig

## Antibody description

Rabbit polyclonal antibody to SIGLEC10

## Preparation

This antigen of this antibody was klh conjugated synthetic peptide derived from human siglec10 321-420/697

## Formulation

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

## Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB

## Dilutions

WB:1:500-2000

## Validations




Fig1: Sample:; HeLa(Human) Cell Lysate at 30 ug; Primary: Anti- SIGLEC10 at 1/1000 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 77 kD; Observed band size: 102 kD