

Anti- antibody

Catalog Number: 176686

Product name

Rabbit IgG

Anti- antibody

Applications

Antibody description

Rabbit polyclonal antibody to

WB, ELISA

Preparation

This antigen of this antibody was recombinant protein


Dilutions

WB: 1:2,000-1:5,000 ELISA 1:1,000- 1:5,000

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Validations



Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Fig1: Western blot analysis on full length recombinant protein of shiga toxin subunit B using polyclonal antibody.

Clonality

Polyclonal

Ig Type