

**Product name**

Rabbit IgG

Anti- antibody

**Applications****Antibody description**

WB, ELISA

Rabbit polyclonal antibody to


**Dilutions****Preparation**

WB: 1:2,000-1:5,000 ELISA 1:1,000- 1:5,000

This antigen of this antibody was recombinant protein

**Validations****Formulation**

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

**Storage**

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

**Clonality**

Polyclonal

**Ig Type**

Fig1: Western blot analysis on full length recombinant protein of shiga toxin subunit B using polyclonal antibody.