

USP48 antibody

Catalog Number: 116686

Product name

USP48 antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

USP48 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in PC-3 cells, HeLa cells. Positive IHC detected in human pancreas cancer tissue. Observed molecular weight by Western-blot: 119 kDa

Preparation

This antibody was obtained by immunization of USP48 recombinant protein (Accession Number: BC013567). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:116686(USP48 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:116686(USP48 antibody) at dilution of 1:100



(under 10x lens)

IP Result of anti-USP48 (IP:Catalog No:116686,

USP48 antibody



Catalog Number: 116686

4ug; Detection:Catalog No:116686 1:300) with HeLa cells lysate 3200ug.