

ABCB7 antibody



Catalog Number: 107686

Product name

ABCB7 antibody

Immunogen

[Human ABCB7 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ABCB7 Rabbit Polyclonal antibody. Positive IP detected in mouse liver tissue. Positive WB detected in mouse liver tissue, mouse thymus tissue. Observed molecular weight by Western-blot: 83kd

Preparation

This antibody was obtained by immunization of ABCB7 recombinant protein (Accession Number: NM_004299). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107686(ABCB7 antibody) at dilution of 1:1000



IP Result of anti-ABCB7 (IP:Catalog No:107686, 4ug; Detection:Catalog No:107686 1:500) with mouse liver tissue lysate 4000ug.