

CACNA2D2 Polyclonal Antibody



Catalog Number: 160386

Product name

CACNA2D2 Polyclonal Antibody

Specificity

Mouse, Rat

Antibody description

Polyclonal antibody to CACNA2D2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 20-200 of human CACNA2D2 (NP_001167522.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:200 - 1:2000

Validations



Western blot - CACNA2D2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using CACNA2D2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 1s.