

# NUP98-NUP96 antibody



Catalog Number: 113386

## Product name

NUP98-NUP96 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NUP98-NUP96 Rabbit Polyclonal antibody.  
Positive WB detected in COLO 320 cells, HeLa cells, Jurkat cells, K-562 cells, MCF7 cells.  
Positive IP detected in COLO 320 cells. Observed molecular weight by Western-blot: 98kd

## Preparation

This antibody was obtained by immunization of NUP98-NUP96 recombinant protein (Accession Number: BC012906). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:113386(NUP98-NUP96 antibody) at dilution of



1:1000

IP Result of anti-NUP98-NUP96 (IP:Catalog No:113386, 4ug; Detection:Catalog No:113386 1:1000) with COLO 320 cells lysate 2000ug.