

Anti-ZNF337 antibody



Catalog Number: 176386

Product name

Anti-ZNF337 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to ZNF337

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human znf337. aa range:471-520

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IF

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

IF: 1:100-1:500

Validations



Fig1.; Western Blot analysis of various cells using ZNF337 Polyclonal Antibody. Secondary antibody [catalog#: HA1001) was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).