

# CLEC4M antibody

Catalog Number: 109386

## Product name

CLEC4M antibody

## Specificity

Human; other species not tested.

## Antibody description

CLEC4M Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue. Positive WB detected in Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of CLEC4M recombinant protein (Accession Number: NM\_214675). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109386(CLEC4M Antibody) at dilution of 1:3000



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109386(CLEC4M Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109386(CLEC4M Antibody) at dilution of 1:50 (under 40x lens)