

# VPS37B antibody

Catalog Number: 116786

---

## Product name

VPS37B antibody

## Specificity

Human; other species not tested.

## Antibody description

VPS37B Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, COLO 320 cells. Observed molecular weight by Western-blot: 31 kDa

## Preparation

This antibody was obtained by immunization of VPS37B recombinant protein (Accession Number: NM\_024667). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116786(VPS37B Antibody) at dilution of 1:1000