

ADH1B antibody



Catalog Number: 107786

Product name

ADH1B antibody

Immunogen

[Human ADH1B Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ADH1B Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 40 kDa

Preparation

This antibody was obtained by immunization of ADH1B recombinant protein (Accession Number: NM_000668). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107786(ADH1B Antibody) at dilution of 1:1000