

Secretogranin III antibody



Catalog Number: 115096

Product name

Secretogranin III antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Secretogranin III Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 50-53 kDa

Preparation

This antibody was obtained by immunization of Secretogranin III recombinant protein (Accession Number: NM_013243). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115096(SCG3 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:115096(SCG3 antibody) at dilution of 1:50 (under 10x lens)